



TS7000 U14 OR Table

UNIVERSAL OPERATING TABLE
FOR GREATER PATIENT ACCESS



TS7000 U14 Mobile Operating Table

Make your operating room more flexible

- All the advantages of the TS7000 OR Table but with a shortened back section, allowing greater access for the anesthesia provider regardless of the surgery.
- The modular design along with a comprehensive line of accessories and components provides the ability to configure the table for a wide range of applications and offers surgeons greater visibility in applications such as gynecology, urology and orthopedics.
- Simplified configuration and setup with improved coupling point system.
- Radiolucent carbon fiber components are also available for full patient intraoperative imaging.



Components & Accessories for Various Applications



Superior access during urology and gynecology procedures with stirrups and URO cut-out.



Simplified patient positioning in urology and gynecology procedures with four-part leg section.



Advanced orthopedic and trauma applications supported by Extension Unit with MIS hip device.



Spine and neuro procedures on universal OR Table with modular structure and broad range of accessories.

ORDERING INFORMATION

	TS7000 U14	
Material-N°	TS7000 U14 (MBW) ¹ Pad TS7000 U14 B	Mat.-No. 2065385 Mat.-No. 2065417
Height Adjustment	59 cm – 114 cm / ± 1 cm (23" – 45" / ± 0.4")	
Trendelenburg/Anti-Trendelenburg	± 45°	
Tilt	± 30°	
Longitudinal Slide	46 cm ± 1 cm (18" / ± 0.4")	
Joint Leg Section Up/Down	+ 90° / - 105°	
Joint Back Section Up/Down	± 90°	
Maximum Load	450 kg (1000 lbs.)	

¹Motorized Base with Wireless Interface

For more information please contact your local Hillrom representative.

hillrom.com

Hillrom
Amerika Building (7th Floor), Hoogoorddreef 15
1101 BA Amsterdam, The Netherlands

TS7000 U14 naming corresponds to the following CE registered products: TruSystem 7000 U (MB) - 1604786; TruSystem 7000 U - 1604788; TruSystem 7000 U (dV) - 1723633; TruSystem 7000 - 1841046; TruSystem 7000 (MBW) - 1841048; TruSystem 7000 (dV) - 1841049; TruSystem 7000 V - 1841050; TruSystem 7000 (MBW) V - 1841082; TruSystem 7000 (dV) V - 1841083; TruSystem 7000 U14 (MBW) - 2065385; TruSystem 7000 U14 (MBW) V - 2065386; Cladding Protection TS7000 1604786 - 2069528; Tam ISI - 1832238; Cable connection TS7000 (dV) - 1816914

The TS7000 U14 Surgical Table is intended to be used by clinicians and medically qualified personnel. This medical device is a regulated healthcare product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. TRUMPF Medizin Systeme GmbH + Co. KG is a subsidiary of Hill-Rom Holdings, Inc.